

Instruction Manual

VDPro® Rabies FA Reagent

CAT.NO. RC-RAB-11



! New version Notification: Read instruction manual before use.

1. Introduction

The **VDPro® Rabies FA reagent** is indirect fluorescent antibody (FA) test reagent. The FA reagent contains pre-diluted monoclonal antibody against nucleocapsid (NC) protein of rabies virus and FITC conjugated anti-mouse conjugate. This reagent is intended for the detection of rabies NC proteins in culture and in acetone-fixed brain and submaxillary tissues of infected animals.

2. Principle

The **VDPro® Rabies FA reagent** uses a pre-diluted monoclonal antibody against the rabies nucleocapsid protein to detect the virus in infected tissue. The antibody is incubated with rabies-infected tissue and will bind to rabies antigen present. Unbound antibody is removed by washing and the antigen-antibody complex is visualized by following reaction of FITC anti-mouse conjugate. Rabies protein in infected cells will fluorescence bright apple-green, background will be stained bright red cause of Evans blue counterstained.

3. Reagents Contents

Anti-Rabies monoclonal antibody(Mab)	-----	8mℓ
Rabies Mab clone 2C6	above 0.1ug/ml	
FITC Anti-mouse Conjugate	-----	8mℓ
FITC Goat Anti-mouse IgG	1/100	
Evans blue counter stain	less than 0.01%	
10X Washing Buffer	-----	120mℓ
FA Mounting Fluid	-----	3mℓ

4. Materials and instruments required (not supplied)

- 1) Slide glass and cover slip
- 2) Cryocut Microtome
- 3) Fluorescent microscope (FITC filter: excitation peak = 490nm, emission peak = 515nm)
- 4) Humid chamber
- 5) Coplin staining jars
- 6) Homogenizer (aerosol-tight)
- 7) Normal and infected mouse brain suspensions
- 8) Positive and negative control slides

5. REAGENTS PREPARATION

1X Washing Buffer

Mix 90mℓ of distilled water with 10mℓ of 10X Washing Buffer.
Keep in refrigerate and use within 2 weeks

PRECUTION

- 1) For *in vitro* use only.
- 2) Do not use past expiry date.
- 3) Incubate the reagent in room temperature before use and return to refrigerator after use.
- 4) Do not allow the reagent to dry on the slides during the staining procedure.
- 5) Handle all specimens, slides and materials coming in contact with them as potentially infectious. Decontaminate with 0.05% sodium hypochlorite.
- 6) Pooling or alteration of any reagent may cause erroneous results.
- 7) Acetone is extremely flammable and harmful if swallowed or inhaled. Keep away from heat, sparks, or flames. Use adequate ventilation and avoid breathing vapor.
- 8) Performance of the fluorescence microscope is critical in achieving satisfactory test results. Microscope objectives, bulb intensity and wattage, and filters may affect results.
- 9) Reagent that is visibly cloudy should not be used.

6. Sample preparation

- ❖ Fluorescent antibody testing is most often performed on brain (hippocampus Ammon's horn) and submaxillary gland specimens. Each rabies specimen should be collected so as to avoid cross-contamination of specimens and contamination of laboratory surfaces.
- ❖ Specimens to be tested should be stored at 2-8°C if testing will be performed within 24 hours. If extended storage is required, specimens should be stored at -70°C or colder.

7. Immunofluorescence Procedure

- 1) Remove the acetone-fixed Positive and Negative Control slides and tissue section slides from the freezer and allow equilibrating to room temperature.
- 2) Add 4-5 drops of 1X anti-rabies monoclonal antibody to slide. Care must be taken to avoid cross contamination of slides during the staining procedure.
- 3) Incubate 30min at room temperature (or 37°C incubator) in humid chamber.
- 4) Remove excess solution from the slides by briefly rinsing with 1X Washing Buffer. Next wash the slide 2-3 times briefly in 1X Washing Buffer.
- 5) Add 4-5 drops of 1X FITC anti-mouse conjugate and incubate 30min at room temperature (or 37°C incubator) in humid chamber.
- 6) Remove excess solution from the slides by briefly rinsing with 1X Washing solution. Next wash the slide 3 times briefly in 1X Washing Buffer.
- 7) Mount cover slip on slides using 2-3 drop of FA mounting fluid.
- 8) Examine the slide on Fluorescent microscope with filter adapted in FITC.

8. Interpretation

- ❖ Specific green fluorescent like large and small dust like inclusions indicate positive reaction.
- ❖ If non-specific reaction was shown, retest using 1:2 diluted FITC anti-mouse conjugate with PBS.
- ❖ Demonstration of non-specific staining (such as that due to certain bacteria) can be confirmed by staining with a FITC-IgG conjugate directed against an unrelated virus, such as canine distemper virus.

Note:

- ❖ Non-specific staining may occur due to non-specific binding of the reagent to leukocytes and certain types of connective tissue. These reactions are morphologically distinct from the reaction seen with rabies-infected cells.
- ❖ The intensity of fluorescence seen in any positive specimen will be a function of the microscope used including filter set and light source, dilution of the Reagent and the quality of the tissue specimen.
- ❖ This reagent will detect rabies-related lyssaviruses other than rabies. Further testing may be necessary to confirm virus type.
- ❖ Slides mounted with glycerol greater than 10% and pH less than 8.0% may be subject to fading and loss of staining intensity.

Median Diagnostics Inc.

878, Sunhwan-daero, Dongnae-myeon, Chuncheon-si,
Gangwon-do, 24399, Republic of Korea
Tel: +82 (0)33 244 0100
Fax: +82 (0)33 244 4634
E-mail: median@mediandx.com