

VDx[®] ASFV qPCR

African Swine Fever Virus real time PCR kit

Advantages of VDx[®] ASFV qPCR

- Confirmatory diagnosis of ASFV
- High sensitivity and specificity
- High reproducibility and high repeatability
- Applicable to all suspect specimens (blood, tissue, feces, livestock processed goods, and food wastes)
- Most conserved p72 gene targeting for detecting all genotypes of ASFV



Performances of VDX[®] ASFV qPCR


- Detection limit: **1 copy/μl**
- Sensitivity: **100%**
- Specificity: **100%**
- Cross reactivity: **0%**

Cut-off

- Less than
Ct-value 40
: **Positive**

testing 19 kinds of bacteria
and viruses including CSFV

Detection limit of VDX[®] ASFV qPCR



Detection limit
: 1 copy/μl

○ Sensitivity for standard DNA of ASFV genotypes

No	Sample (1/10x)	Country of origin	Year	P72 genotype	VDX® qPCR		No	Sample (1/10x)	Country of origin	Year	P72 genotype	VDX® qPCR	
					Ct	Result						Ct	Result
1	E70	Spain	1970	I	27.69	POS	13	MOL16/DP-CERNO1	Moldova	2016	II	31.7	POS
2	BF07/IqTC	Burkina Faso	2007	I	25.71	POS	14	MOL16/DP-MOSA1	Moldova	2016	II	32.01	POS
3	SS14/WB-Sassari1	Italy	2014	I	32.72	POS	15	Moz64	Mozambique	1964	V	27.37	POS
4	SS14/DP-Cagliari1	Italy	2014	I	35.44	POS	16	MwLil 20/1	Malawi	1983	VIII	31.7	POS
5	Arm07	Armenia	2007	II	28.99	POS	17	Ken11/KisP52	Kenya	2011	IX	30.41	POS
6	Ukr12/Zapo	Ukraine	2012	II	29.05	POS	18	Ken06.Bus	Kenya	2006	IX	29.54	POS
7	Ukr15/DP-Kieve 1	Ukraine	2015	II	29.24	POS	19	Ken08Tk.2/1	Kenya	2007	X	28.27	POS
8	LT14/1490	Lithuania	2014	II	29.88	POS	20	UG10/Tk3.2	Uganda	2010	X	33.81	POS
9	Pol14/Krus	Poland	2014	II	29.74	POS	21	Eth13/1505	Ethiopia	2013	XXIII	31.24	POS
10	Lv14/DP/Robez3	Latvia	2014	II	32.89	POS	22	PC				22.4	
11	Est14/WB-Valga-1	Estonia	2014	II	32.5	POS	23	NC				N/A	
12	Est15/Wb-Tartu14	Estonia	2015	II	32.23	POS							

○ Sensitivity and specificity for ASFV-positive and negative tissues and sera

No	Sample	Clinical	Isolate	Genotype	OIE qPCR	VDX® qPCR(Ct)	No	Sample	Clinical	Isolate	Genotype	OIE qPCR	VDX® qPCR(Ct)
1	S17	Acute	Ukr12/Zapo	II	3+	25.03	11	S27	Tonsil			-	-
2	S18	Chronic	NH/P68	I	+	33.74	12	S28	Acute	Ukr12/Zapo	II	3+	26.24
3	S19	Acute	L60	I	3+	24.2	13	S29	S28, 1/50x			3+	24.55
4	S20	Subacute	Ken05/Tk1	X	3+	23.52	14	S30	S28, 1/400x			2+	28.21
5	S21	Acute	Arm07	II	3+	24.48	15	S31(Serum)	Serum			-	-
6	S22	Acute	Ken06.bus	IX	3+	23.11	16	S32	Acute	Ukr12/Zapo	II	+	29.07
7	S23	S17, 1/200x			+	33.16	17	PC					22.97
8	S24	S22, 1/200x			2+	31.06	18	NC					-
9	S25	Kidney			-	-	19						
10	S26	Lung			-	-	20						

Cross reactivity of VDX® ASFV qPCR with 19 different pathogens (12 viruses, 7 bacteria)

No	Pathogens	Strain	Concentration	Result
1	PRRSV	PL96-1	10 ⁴ TCID ₅₀ /ml	NEG
2	PRRSV	Lelystad	10 ² TCID ₅₀ /ml	NEG
3	EMCV	VR129	10 ^{6.5} TCID ₅₀ /ml	NEG
4	JEV	Nakayama	10 ^{4.8} TCID ₅₀ /ml	NEG
5	CSFV	LOM	10 ⁵ TCID ₅₀ /ml	NEG
6	SIV	H1N2	10 ⁵ TCID ₅₀ /ml	NEG
7	ADV	YS	10 ⁵ TCID ₅₀ /ml	NEG
8	PPV	PV9	10 ⁵ TCID ₅₀ /ml	NEG
9	PCV2	field isolation	10 ⁵ TCID ₅₀ /ml	NEG
10	PEDV	SM98	10 ⁵ TCOD ₅₀ /ml	NEG
11	TGEV	Pyeong Taek	10 ⁴ TCOD ₅₀ /ml	NEG
12	Rotavirus	OSU	2x10 ³ TCOD ₅₀ /ml	NEG
13	App 2	ATCC27089	3.2X10 ⁶ CFU/ml	NEG
14	App 5	ATCC33377	2X10 ⁶ CFU/ml	NEG
15	P.m._A	ATCC43137	2.1X10 ⁶ CFU/ml	NEG
16	P.m._D	ATCC12948	2.2X10 ⁶ CFU/ml	NEG
17	M.hyo	J	-	NEG
18	H.p.	ATCC19417	4.28X10 ⁵ CFU/ml	NEG
19	Sal. typhi.	Sal13	-	NEG