



PRODUCT SPECIFICATION

CATALOG NUMBER : 22304

PRODUCT	rec. Malaria Pv MSP 1-2 [<i>Plasmodium vivax</i> merozoite surface protein]
PRODUCT DESCRIPTION	The recombinant Pv-MSP 1-2 was expressed in <i>E. coli</i> and purified by affinity chromatography to homogeneity. The purified recombinant protein was concentrated and dialyzed with phosphate-buffered saline.
CATALOG No.	22304
EXPRESSED IN	<i>E. coli</i>
PRESERVATIVES	0.05% Sodium Azide
BUFFER	10mM PBS, pH 7.2
BIOHAZARD	Use normal laboratory precautions and procedures.
APPLICATIONS	Rapid Diagnostic Test(RDT), EIA
PURITY	>95% by SDS-PAGE gel
STORAGE	Store at -20°C.

For In Vitro Use Only