

## **POLYCLONAL ANTIBODY SPECIFICATION**

**CAT. No. 9904**

**Products : Anti-PCV2 Rabbit Poly-Antibody**

**Descriptions:** Purified rabbit antibody immunized with PCV2 recombinant NC protein.

**Specificity:** Reactive to Porcine Circovirus Type 2.  
Low level cross reactive to PCV Type 1 in Fluorescent antibody(FA) test.

**Immunogen:** Purified recombinant NC protein of PCV2 (rNC).

**Host animal:** Serum from Rabbit, Newzealand white.

**Purification:** Protein A affinity chromatography.

**Buffer:** PBS (pH 7.4) with 0.1% Sodium azide.

**Applications:** Fluorescent Antibody (FA) Test,  
Enzyme Linked Immnosorbent Assay (ELISA).  
Optimum working dilutions must be determined by user.

**Storage:** Store at 2~8°C,  
Antibody solutions can be stored at -20°C but avoid repeated freezing and thawing.

*In Vitro use only.*  
*Research use only*